


High Quality Custom

Antisense - RNAi

Oligonucleotides


A strong parallel approaches for your gene silencing expedition.



Why should you accept less? Choose from a complete range of quality ready-2-use custom RNA and duplex RNAi services delivered, when and where they're needed from a trusted brand with fully deprotected, desalted and MS controlled products.

Looking for a laboratory to manufacture your RNA with PTOs for antisense experiments, or with a variety of other modifications and multiply your options for gene silencing? visit us at

www.biosyn.com

- 
-  *High throughput synthesis*
 -  *Large scale productions*
 -  *Various highly specific modifications, service options*
 -  *Dual fluorophore labelings*
 -  *ASR or cGMP oligonucleotide manufacturing*



The most trusted name in Quality Custom Synthesis



THE HOME FOR RNA SCIENTISTS

The **RNA Society** is an interdisciplinary, cohesive intellectual home for those interested in all aspects of RNA Science.

We welcome new members from all areas of scientific research and we look forward to sharing the new perspectives they bring to the Society. There are many benefits to being a member of the **RNA Society**.

These include:

- Print copies and web access to the Society journal, **RNA** (Impact factor of 6.1 and an excellent immediacy index)
- Reduced costs for publishing your papers in **RNA**.
- The RNA Society Newsletter, a forum for disseminating information to members and discussing issues affecting the Society and RNA Science
- Reduced registration fees for the annual meeting of the Society
- Numerous opportunities for junior scientists to become involved in the Society
- The Directory of Members, available on the web and in print
- Free job postings on the Society Employment and Careers website
- Opportunities to request Travel Fellowship and Meeting Support for RNA-related meetings you are organizing.

**VISIT US AT WWW.RNASOCIETY.ORG AND BECOME
A MEMBER OF THE RNA SOCIETY TODAY!**

THE RNA SOCIETY, 9650 ROCKVILLE PIKE, BETHESDA, MD 20814-3998 USA

PHONE: 301-634-7120 FAX: 301-634-7420

EMAIL: RNA@FASEB.ORG

[HTTP://WWW.RNASOCIETY.ORG/](http://WWW.RNASOCIETY.ORG/)



siRNA delivery to neuronal cells. What a wonderful world.

Thermo Scientific Dharmacon Accell siRNA provides a breakthrough in siRNA delivery. It is specially modified for delivery without a separate transfection reagent, and has demonstrated excellent results for gene silencing in numerous neuronal cell types. RNAi in neuroscience research is no longer constrained by conventional delivery methods that rely on lipid-based reagents, viruses, or instrumentation.

- Ex vivo mouse brain slice cultures
- Primary rat cortical neurons
- Human neural stem cells

Learn More Today

For more information, call 800-235-9880 or check our website for additional information:
www.thermo.com/acell

Visit us at the
Society of
Neuroscience
Booth #2301

**Dharmacon® Accell™ siRNA can
accelerate discovery**

Stop by our booth to learn more about
this breakthrough RNAi technology.

Moving science forward

Thermo
S C I E N T I F I C

Part of Thermo Fisher Scientific

Save
15% in
October!

Ina K. Dahlsveen, Ph.D.

Start your microRNA discovery from a single cell


With the new miRCURY LNA™ microRNA PCR System, your discovery starts from 10 pg total RNA – you only need a single cell. You can even detect microRNAs with high and low expression from the same sample in your experiment, as the dynamic range is 8 logs wide.

Our PCR system combines microRNA-specific primers and SYBR® Green detection and is part of Exiqon's popular miRCURY LNA™ product line, the most complete range of tools for microRNA analysis available today.

Ina K. Dahlsveen, Ph.D.,
Product Manager



Accurate quantitation over a wide dynamic range of microRNA copy numbers



Amplification of a synthetic microRNA was performed using the miRCURY LNA™ microRNA PCR System. Hsa-miR-100 (2×10^9 copies) input was used to generate cDNA. A ten-fold serial dilution ranging from 2×10^9 copies to 2 copies cDNA was amplified by real-time PCR. **A**, the amplification curves and **B**, the standard curve show excellent linear correlation between the decreasing cycle numbers and the logarithm of the microRNA copy number, proving the assay enables sensitive microRNA detection down to less than 10 copies with a wide dynamic range of at least 8 logs.

Get the complete miRCURY LNA™ System now
and save 15%! Visit exiqon.com/PCR15

Offer only valid in October 2008

EXIQON
Seek Find Verify