

Study m6A Modifications Like Never Before




Fig. 1. m6A modified RNAs are immunoprecipitated and compared with the input RNAs in two-color channels of the Epitranscriptomic Array to accurately measure the modification percentage and abundance level for each transcript. The percentage of modification of “transcript a” RNA determines its different cell fates, e.g. for protein translation or RNA decay.

Epitranscriptomic Array

- Quantifying the percentage of modification for each transcript
- Coverage of mRNAs, lncRNAs, circRNAs and more
- RNA modifications at transcript-specific level
- Low sample amount, starting from 1 ug total RNA




Fig. 2. Unmethylated (ACA) sites, but not methylated (m6ACA) sites, are cleaved by RNase MazF. The MazF treated and corresponding untreated sites are compared in two-color channels of the m6A Single Nucleotide Array to accurately quantify the m6A modification percentage, and abundance level at precise single nucleotide resolution.

m6A Single Nucleotide Array

- A new orthogonal methodology for m6A detection
- Single-Nucleotide resolution for m6A site location
- Quantifying the percentage of m6A site modification
- Reliable m6A site collection and systematic annotations
- Low sample amount, starting from 1 ug total RNA

Increased Viral RNA Yields for the Earliest Detection of COVID-19 Outbreaks

Zymo Environ Water RNA Kit For SARS-CoV-2 Wastewater Surveillance



Highest RNA Concentration
Removes PCR Inhibitors
Eliminates Infection Risk

The Zymo Environ™ Water RNA Kit simplifies COVID-19 wastewater surveillance by increasing viral enrichment, inactivating pathogens, and improving limits of detection. The workflow yields highly concentrated, PCR inhibitor-free RNA ready for downstream applications including RT-qPCR, RT-ddPCR, and NGS.

Upgrade Your COVID-19 Wastewater Surveillance Today!

Visit www.zymoresearch.com/pages/zymo-environ-water-rna-kit



www.zymoresearch.com



info@zymoresearch.com



Toll Free: (888) 882-9682



Capture the entire transcriptome— in a single library.

Introducing the **new** SEQuoia Complete Stranded RNA Library Prep Kit for NGS.



It's new. It's different. And it's only from Bio-Rad.

The SEQuoia Complete Stranded RNA Library Prep kit gives you a truly holistic view of the transcriptome, capturing all types of RNAs (long, short, and everything in between) from all types of samples in less than 4 hours. With unparalleled uniformity of coverage, strandedness, and efficiency, what you will discover will definitely surprise you.

Are you ready for a new approach to transcriptomics?
Get the details at [bio-rad.com/SEQuoiaCompleteRNA](https://www.bio-rad.com/SEQuoiaCompleteRNA)

BIO-RAD

#ScienceForward

AACR

American Association
for Cancer Research®

ANNUAL MEETING

2021

VIRTUAL ★ #AACR21



APRIL 10-15
MAY 17-21

**We look forward to
seeing you in 2021!**

Register now for the world's most comprehensive
annual meeting dedicated to the research,
prevention, detection, and treatment of cancer.

★ **Become a Member!**

Join the AACR and receive a discount on registration.

**VISIT AACR.ORG/AACR2021 FOR MORE
INFORMATION AND TO REGISTER!**

AACR.ORG ★ #AACR21



AACR.ORG/AACR2021



Continuing Medical Education Activity -
AMA PRA Category 1 Credits™ available


DISCOVERY
SCIENCE
DRIVING
CLINICAL
BREAKTHROUGHS

Meet your custom RNA requirements with our nearly limitless synthetic options.

If your NMR, crystallography, or RNA binding study requires a non-standard chemical modification or a commercially unavailable modified nucleobase, we can help. For more than 20 years our chemists have utilized Dharmacon™ proprietary 2' ACE chemistry to synthesize RNA with superior yield and quality giving you the flexibility you need.

Better RNA yields give you more for your money.

RNA yield (HPLC-purified)



The yield for RNA oligos of three different lengths was compared for all available synthesis scales between Dharmacon and a competitor.