



Functionality,
reimagined.



The New Vortex Mixer

- Mix and vortex in one compact device
- Precisely set speed and duration controls
- Glove compatible touchscreen
- Quickly vortex tubes with just one touch

The Vortex Mixer is ideally suited for assay protocols that require gentle mixing, all the way up to resuspensions that require vigorous vortexing. Find out more at usascientific.com/vortex_mixer.



Schedule a
demo!

USn
SCIENTIFIC
SIDE BY SCIENCE



Non-coding RNA and Epitranscriptomic Solutions



Get Your Free eBook

Circular RNA Arrays

Accurately profile circular RNAs by highly specific circular junction probe design

LncRNA Arrays

Overcome the limitations of RNA-seq for lncRNAs often at low abundance

Small RNA Arrays

Accurately profile miRNA, pre-miRNA, tRNA, tsRNA, and snoRNA simultaneously

Epitranscriptomic Arrays

Quantify the percentage of m6A modifications at the transcript specific level


m6A Single Nucleotide Arrays

Locate and quantify the exact m6A site at single nucleotide resolution

Pioneering Excellence in Immunotherapy Research

ProMab Biotechnologies' offers a one-stop-shop for contract research services to enhance antibody, mRNA-LNP, and cell therapy research.

As the world's first CAR-T CRO and a long-established antibody platform for over 20 years, ProMab can leverage diverse innovation pathways to fuel your next therapeutic research adventure.



mRNA-Lipid Nanoparticles

CAR-T/NK/Macrophage

Monoclonal & Bispecific Antibodies


ProMab Biotechnologies' combines advanced synergies under one umbrella:

- Completely validated mRNA-LNP platform supporting formulation and manufacturing
- Chimeric antigen receptor design and development for a plethora of immune cell types
- Flexible antibody production and engineering offerings for research and therapeutic applications

Imagine the possibilities partnering with our end-to-end services.
For more details, scan the QR code or visit www.promab.com




2600 Hilltop Dr, Building B, Richmond, CA 94806
1.866.339.0871 | info@promab.com



AACR American Association
for Cancer Research®

2024 SCIENTIFIC CONFERENCES

Presenting the most significant research on cancer etiology, prevention, diagnosis, and treatment



Bladder Cancer: Transforming the Field May 17-20, 2024 | Charlotte, NC

Conference Cochair: Lars Dyrskjøjt Andersen,
Donna E. Hansel, Dan Theodorescu, and
Tahlita C. M. Zuiverloon

Expanding and Translating Cancer Synthetic Vulnerabilities June 10-13, 2024 | Montreal, Quebec, Canada

Conference Cochair: Kimberly Stegmaier,
E. Alejandro Sweet-Cordero, Michael Erb, and
Kris C. Wood

Fourth AACR International Meeting Advances in Malignant Lymphoma: Maximizing the Basic-Translational Interface for Clinical Application In Cooperation with the International Conference on Malignant Lymphoma (ICML) In Association with *Blood Cancer Discovery* June 19 - 22, 2024 | Philadelphia, PA

Conference Chair: Michael R. Green
Conference Cochair: Francesco Bertoni*,
Christopher R. Flowers, Ari M. Melnick,
Laura Pasqualucci, and Margaret A. Shipp
*International Conference on Malignant Lymphoma (ICML)
Representative

Advances in Pediatric Cancer Research In Collaboration with the AACR Pediatric Cancer Working Group (PCWG) September 5-8, 2024 | Toronto, ON, Canada

Conference Cochair: Alejandro Gutierrez,
Cynthia E. Hawkins, Andrea A. Hayes, and
Gilles Vassal

Pancreatic Cancer September 15-18, 2024 | Boston, MA

Conference Cochair: Peter J. Allen, Stephanie K.
Dogan, Michael A. (Tony) Hollingsworth,
and Alec C. Kimmelman

17th AACR Conference on the Science of Cancer Health Disparities in Racial/Ethnic Minorities and the Medically Underserved In Association with the AACR Minorities in Cancer Research Council

September 21 - 24, 2024 | Los Angeles, CA
Conference Cochair: Moon S. Chen, Melissa
B. Davis, Carmen E. Guerra, Ruben A. Mesa,
Jasmine Plummer and Susan M. Shinagawa

Tumor Immunology and Immunotherapy In Association with the AACR Cancer Immunology Working Group (CIMM) October 18-21, 2024 | Boston, MA

Conference Cochair: Yvonne Y. Chen, Sergio
Quezada, Robert D. Schreiber, and Fernando
Vidal-Vanaclocha

EORTC-NCI-AACR Symposium on Molecular Targets and Cancer Therapeutics October 23-25, 2024 | Barcelona, Spain

Conference Cochair: E.G. Elisabeth de Vries,
Tim F. Greten, and Timothy A. Yap

Liquid Biopsy: From Discovery to Clinical Implementation November 13-16, 2024 | San Diego, CA

Conference Cochair: Catherine Alix-Panabières,
Luis A. Diaz, Jr., Maximilian Diehn, Y.M. Dennis Lo,
and Klaus Pantel

RNAs as Drivers, Targets, and Therapeutics in Cancer November 14-17, 2024 | Bellevue, WA

Conference Cochair: Howard Y. Chang, Joshua T.
Mendell, Anastasia Khvorova, and V. Narry Kim

Learn more and register at
AACR.org/Calendar


Follow us @AACR



2401078C

GENE TOOLS, LLC

www.gene-tools.com



A Morpholino oligo works like sequence-specific masking tape, blocking a section of RNA.

Halt the initiation complex.

Block snRNP binding.

Inhibit miRNA maturation.

Mask a poly-A tail signal.

Switch splice regulation.

Invade RNA secondary structure.

Proven technology:

Four FDA-approved drugs

Over 11,000 publications

Nurture Growth, Harvest Discoveries

Unveil the Complete Plant
Transcriptome with RiboCop

Enhancing Plant Transcriptomics

Get the best sequencing results with convenient end-to-end RNA-Seq workflows for plant transcriptomics, expression profiling and plant-microbe interaction studies!

Key Benefits for Plant RNA-Seq with RiboCop



Unlock the full potential of low input and lowly concentrated samples - RiboCop supports a wide input range from **1 ng - 1 µg** total RNA.



Enzyme-free rRNA depletion preserves full-length RNA and is ideally suited for processing even degraded samples.



Simultaneous depletion of plant and microbe with convenient combinations of RiboCop for Plants and RiboCop Bacteria kits.

Plant total RNA is comprised of large amounts of undesired ribosomal RNA (rRNA) which accounts for up to ~97 % of all transcripts. Lexogen's RiboCop rRNA Depletion Kits for Plants **remove cytoplasmic (5.8S, 18S and 25S), mitochondrial (5S, 18S and 26S) and plastid (4.5S, 5S, 16S and 23S) rRNA** from intact as well as degraded input RNA. Ribodepletion provides an unbiased view of the transcriptome and allows simultaneous studies of plant and microbial RNA profiles. RiboCop for Plants has been validated for various species, including *A. thaliana*, *Z. mays*, *G. max*, and *H. lumulus*. RiboCop is available for various species - find your match today!

RiboCop rRNA Depletion Kits - Find Your Match!

Cat. №	Product Name
rRNA Depletion for Plants	
237	RiboCop rRNA Depletion Kits for Plants
rRNA and Globin Depletion for Human/Mouse/Rat	
144	RiboCop rRNA Depletion Kits for Human/Mouse/Rat
145	RiboCop rRNA Depletion Kits for Human/Mouse/Rat Plus Globin
rRNA Depletion for Fish	
241	RiboCop rRNA Depletion Kits for Fish
rRNA Depletion for Bacteria	
125	RiboCop rRNA Depletion Kits for Mixed Bacterial Samples (META)
126	RiboCop rRNA Depletion Kits for Gram Negative Bacteria (G-)
127	RiboCop rRNA Depletion Kits for Gram Positive Bacteria (G+)
rRNA Depletion for Yeast and Fungi	
190	RiboCop rRNA Depletion Kits for Yeast

Interested to learn more?

Scan the QR-code or visit our [webpage](#) for more information.

