

Supplemental_Figure_S6

dsRNA immunoblot showing reduction of 3'-extended dsRNA byproducts with TsT7-2.

Immunoblot using dsRNA antibody J2 on *in vitro* transcription reactions with CLuc template using TsT7-2 at 37°C or 50°C for 1 hour. *In vitro* transcription reactions using IVT Template short 1 (30 bp) ran on denaturing gels under standard conditions (four rNTPs), and under conditions where one rNTP is omitted to prevent formation of the dsRNA products (-rUTP).

