

Supplemental Information

Desthiobiotin (DTB)-modified and TAMRA-modified 2''OMeNAD⁺ are RNA 2'-phosphotransferase (Tpt1) poisons that enable affinity-tagging and fluorescence-tagging of internal RNA 2'-phosphate groups

Renata Kasprzyk, Shreya Ghosh, and Stewart Shuman

Supplemental Methods

Figure S1

HPLC profiles, HR MS and NMR spectra

General methods for chemical synthesis

All solvents and chemical reagents used for compound synthesis were purchased from Sigma-Aldrich, VWR Chemicals, and Acros Organics and were used without any prior treatment. The synthesis of triethylammonium salt of 1-(2'-methoxy- β -D-ribofuranosyl)nicotinamide-5'-O-monophosphate (2"OMeNMN) has been reported in our previous work (Arnold et al. 2025). Triethylammonium salt of N6-(6-aminohexyl)adenosine 5'-monophosphate (N6-hex-AMP) was synthesized as described previously (Szczepaniak et al. 2012). Desthiobiotin NHS ester (DTB-NHS), *P*-imidazolides of adenosine 5'-monophosphates modified with desthiobiotin and TAMRA at the C2 position (C2-DTB-AMP-Im and C2-TAMRA-AMP-Im), their precursors, NAD-2-DTB, and NAD-2-TAMRA were synthesized as described previously (Kasprzyk et al. 2024; Lehner et al. 2022).

The progress of nucleotide synthesis was monitored using analytical RP-HPLC. The analysis was performed using a Shimadzu system consisting of LC-20AT pumps, an SPD-M20A absorption detector, a SIL-20A autosampler, and a CBM-20A communication module. The system was equipped with an EC HPLC NUCLEODUR C18 Pyramid column (5 μ m, 250 x 4.6 mm, flow rate of 1.0 mL/min), operating with a linear gradient of acetonitrile in 0.1 % formic acid. UV-detection was conducted at 254 nm, with additional detection at 555 nm for TAMRA-labelled nucleotides.

Compounds **1**, **3**, and **4** (Fig. S1) were purified using semipreparative RP-HPLC on a Shimadzu system consisting of LC-20AP pumps, an SPD-M20A absorption detector, an FRC-10A fraction collector, and a CBM-20A communication module. The system was equipped with a VP HPLC NUCLEODUR C18 Pyramid column (5 μ m, 250 x 16 mm, flow rate 8.0 mL/min) using different buffers and eluents (see the table below). UV-detection was performed at 254 nm, with additional detection at 555 nm for TAMRA-labelled nucleotides. The specific HPLC methods used in the studies are outlined below:

Method name	Eluent A	Eluent B	Eluent gradient		Column
P1	0.05 M triethylammonium acetate buffer (pH ~ 6)	ACN	5 %	0 – 10 min	VP HPLC NUCLEODUR C18 HTec, 5 μ m, 250 x 21 mm
			5 – 80 %	5 – 60 min	
			80 – 100 %	60 – 65 min	
			100 %	65 – 70 min	
			100 – 5 %	70 – 75 min	
			5 %	75 – 80 min	
P2	0.1 % formic acid	ACN	0 %	0 – 5 min	VP HPLC NUCLEODUR C18 Pyramid, 5 μ m, 250 x 16 mm
			0 – 55 %	5 – 40 min	
			55 – 100 %	40 – 45 min	
			100 %	45 – 50 min	
			100 – 0 %	50 – 55 min	
			0 %	55 – 60 min	
P3	0.1 % formic acid	ACN	0 %	0 – 5 min	VP HPLC NUCLEODUR C18 Pyramid, 5 μ m, 250 x 16 mm
			0 – 45 %	5 – 40 min	
			45 – 100 %	40 – 45 min	
			100 %	45 – 50 min	
			100 – 0 %	50 – 55 min	
			0 %	55 – 60 min	
P4	0.05 M ammonium acetate buffer (pH 5.9)	ACN	0 %	0 – 5 min	VP HPLC
			0 – 65 %	5 – 60 min	NUCLEODUR C18
			65 – 100 %	60 – 70 min	Pyramid, 5 μ m,

P5	0.05 M ammonium acetate buffer (pH 5.9)	ACN	100 %	70 – 75 min	250 x 16 mm
			100 – 0 %	75 – 80 min	
			0 %	80 – 85 min	
			0 %	0 – 5 min	VP HPLC NUCLEODUR C18 Pyramid, 5 µm, 250 x 16 mm
			0 – 45 %	5 – 65 min	
			45 – 100 %	65 – 70 min	
			100 %	70 – 75 min	
			100 – 0 %	75 – 80 min	
			0 %	80 – 85 min	

Yields after purification were calculated based on sample absorption, measured using the ND-1000 NanoDrop spectrophotometer (PEQLAB). The structure and homogeneity of the final product were confirmed using a combination of analytical techniques, including RP-HPLC, high-resolution mass spectrometry (HRMS) in ESI+ mode, and various NMR methods: ¹H NMR, ¹H-¹H COSY, ¹³C NMR, ¹H-¹³C HSQC NMR, ³¹P NMR, and ¹H-³¹P HMBC NMR spectroscopy. High resolution mass spectra were recorded on an LTQ Orbitrap Discovery. Routine NMR spectra were recorded at 25°C with nucleotide concentrations ranging from 5 to 10 mM, using a Bruker AVANCE NEO spectrometer operating at 800 MHz (¹H, ¹H-¹H COSY), 201 MHz (¹³C), and 324 MHz (³¹P), a Bruker Avance NEO 600 spectrometer operating at 151 MHz (¹³C), a Jeol ECZ 500R spectrometer operating at 500 MHz (¹H, ¹H-¹H COSY), 126 MHz (¹³C), and 202 MHz (³¹P) or a Bruker Avance III 400 spectrometer operating at 400 MHz (¹H, ¹H-¹H COSY). The ¹H chemical shifts were referenced to the solvent signal, and the ³¹P and ¹³C chemical shifts were adjusted based on correlation spectra with the ¹H signals.

References:

- Arnold J, Ghosh S, Kasprzyk R, Brakonier M, Hanna M, Marx A, Shuman S. 2024. Chemical synthesis of 2"OMeNAD⁺ and its deployment as an RNA 2'-phosphotransferase (Tpt1) 'poison' that traps the enzyme on its abortive RNA-2'-PO₄ -(ADP-2"OMe-ribose) reaction intermediate. *Nucleic Acids Res* 23: 10533-10542.
- Kasprzyk R, Rieth S, Heid P, Stengel F, Marx A. 2024. Cell-permeable nicotinamide adenine dinucleotides for exploration of cellular protein ADP-ribosylation. *Angew Chem Int Ed* 63: e202411203.
- Lehner M, Rieth S, Höllmüller E, Spliesgar D, Mertes B, Stengel F, Marx A. 2022. Profiling of the ADP-ribosylome in living cells. *Angew Chem Int Ed* 61: e202200977.
- Szczepaniak SA, Zuberek J, Darzynkiewicz E, Kufel J, Jemielty J. (2012. Affinity resins containing enzymatically resistant mRNA cap analogs—a new tool for the analysis of cap-binding proteins. *RNA* 18: 1421-1432.






Figure S1. Synthesis of desthiobiotin (DTB)-modified and 5-carboxytetramethylrhodamine (TAMRA)-modified 2''OMeNAD⁺ analogs. (A) (i) N-Hydroxysuccinimide ester of DTB (DTB-NHS), 50 mM borate buffer pH 8.2, r.t., 2 h, 75%. (ii) imidazole, 2,2'-dithiodipyridine, triethylamine, triphenylphosphine, dimethylformamide, r.t., 18 h, then NaClO₄, acetone, 4°C, precipitation, 73%. (B) (i) formamide/dimethyl sulfoxide (1:1, v/v), MgCl₂, r.t., 24 h, 35%. (ii) formamide/dimethyl sulfoxide (1:1, v/v), MgCl₂, r.t., 20 h, 42% for 2''OMeNAD-2-DTB or formamide/dimethyl sulfoxide (3:2, v/v), MgCl₂, MnCl₂, r.t., 22h, 42% for 2''OMeNAD-2-TAMRA.


HPLC profiles, HR MS and NMR spectra




¹H NMR (500 MHz, D₂O, 25°C)




¹H-¹H COSY (D₂O, 25°C)





$^{31}\text{PCPD}$ NMR (202 MHz, D_2O , 25°C)



^1H - ^{31}P HMBC (D_2O , 25°C)




^{13}C NMR 125 MHz, D_2O , 25°C)


 ${}^1\text{H}$ - ${}^{13}\text{C}$ HSQC (D_2O , 25°C)

N6-DTB-AMP-Im (2)


Formula




HRMS (+)ES




¹H NMR (800 MHz, MeOD, 25°C)




^1H - ^1H COSY (MeOD, 25°C)




^{31}P CPD NMR (324 MHz, MeOD, 25°C)




^1H - ^{31}P HMBC (MeOD, 25°C)



^{13}C NMR 201 MHz, MeOD, 25°C




^1H - ^{13}C HSQC (MeOD, 25°C)




2''OMeNAD-6-DTB (3)


Formula




HPLC purification profile (Method P2)




HRMS (+)ES




¹H NMR (400 MHz, D₂O, 25°C)¹H-¹H COSY (D₂O, 25°C)


^{31}P CPD NMR (162 MHz, D_2O , 25°C)




^1H - ^{31}P HMBC (D_2O , 25°C)



^{13}C NMR 151 MHz, D_2O , 25°C)




^1H - ^{13}C HSQC (D_2O , 25°C)




2"OMeNAD-2-DTB (4)


Formula



HPLC purification profile (Method P2)




HPLC repurification profile (Method P3)




HRMS (+)ES


250320_RK_2DTB #1 RT: 0.01 AV: 1 NL: 1.21E6
T: FTMS + p ESI Full ms [800.00-1200.00]




¹H NMR (800 MHz, D₂O, 25°C)




^1H - ^1H COSY (D_2O , 25°C)




^{31}P CPD NMR (324 MHz, D_2O , 25°C)




¹H-³¹P HMBC (D₂O, 25°C)






¹³C NMR 151 MHz, D₂O, 25°C




^1H - ^{13}C HSQC (D_2O , 25°C)




2"OMeNAD-2-TAMRA (5)	
Formula	
HPLC purification profile (Method P4)	
	

HRMS (+)ES


250508_2_37_1#2-99 RT: 0.02-1.00 AV: 98 NL: 5.51E6
T: FTMS + p ESI Full ms [550.00-650.00]




¹H NMR (600 MHz, CD₃COOD, 25°C)




¹H-¹H COSY (CD₃COOD, 25°C)




³¹PCPD NMR (324 MHz, CD₃COOD, 25°C)



^1H - ^{31}P HMBC (CD_3COOD , 25°C)



^{13}C NMR 151 MHz, CD_3COOD , 25°C)



^1H - ^{13}C HSQC (CD_3COOD , 25°C)

